


## Seguinocide F

## Chemical Properties

CAS No. :	220778-37-0
Formula:	C <sub>26</sub> H <sub>32</sub> O <sub>15</sub>
Molecular Weight:	584.527
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Seguinocide F is a useful organic compound for research related to life sciences. The catalog number is T125314 and the CAS number is 220778-37-0.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7108 mL	8.5539 mL	17.1078 mL
5 mM	0.3422 mL	1.7108 mL	3.4216 mL
10 mM	0.1711 mL	0.8554 mL	1.7108 mL
50 mM	0.0342 mL	0.1711 mL	0.3422 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481