


(R)-VX-984

Chemical Properties

CAS No. :	2448475-19-0
Formula:	C ₂₃ H ₂₁ D ₂ N ₇ O
Molecular Weight:	415.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(R)-VX-984 is the (R)-enantiomer of VX-984. VX-984 is a potent inhibitor of DNA-PK.
Targets(IC50)	DNA-PK

Solubility Information

Solubility	DMSO: 11.36 mg/mL (27.34 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 1 mg/mL (2.41 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4068 mL	12.034 mL	24.068 mL
5 mM	0.4814 mL	2.4068 mL	4.8136 mL
10 mM	0.2407 mL	1.2034 mL	2.4068 mL
50 mM	0.0481 mL	0.2407 mL	0.4814 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481