


Yadanzioside P

Chemical Properties

CAS No. : 79439-84-2
 Formula: C₃₄H₄₆O₁₆
 Molecular Weight: 710.726
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Yadanzioside P has antileukemic and antitumor activities.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.407 mL	7.035 mL	14.070 mL
5 mM	0.2814 mL	1.407 mL	2.814 mL
10 mM	0.1407 mL	0.7035 mL	1.407 mL
50 mM	0.0281 mL	0.1407 mL	0.2814 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481