

Palbociclib-SMCC


Chemical Properties

CAS No. :

Formula: C36H42N8O5

Molecular Weight: 666.78

Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Palbociclib-SMCC is a Palbociclib derivative with a SMCC linker.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4997 mL	7.4987 mL	14.9975 mL
5 mM	0.2999 mL	1.4997 mL	2.9995 mL
10 mM	0.150 mL	0.7499 mL	1.4997 mL
50 mM	0.030 mL	0.150 mL	0.2999 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481