


Nivasorexant

Chemical Properties

CAS No. : 1435480-40-2
 Formula: C₂₃H₂₃N₇O₂
 Molecular Weight: 429.47
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nivasorexant is a potent antagonist of the orexin receptor.
Targets(IC50)	Others,OX Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3285 mL	11.6423 mL	23.2845 mL
5 mM	0.4657 mL	2.3285 mL	4.6569 mL
10 mM	0.2328 mL	1.1642 mL	2.3285 mL
50 mM	0.0466 mL	0.2328 mL	0.4657 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481