


Dammarenediol II 3-O-caffeate

Chemical Properties

CAS No. :	171438-55-4
Formula:	C ₃₉ H ₅₈ O ₅
Molecular Weight:	606.888
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dammarenediol II 3-O-caffeate is a natural product from <i>Ostryopsis nobilis</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6477 mL	8.2387 mL	16.4775 mL
5 mM	0.3295 mL	1.6477 mL	3.2955 mL
10 mM	0.1648 mL	0.8239 mL	1.6477 mL
50 mM	0.033 mL	0.1648 mL	0.3295 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481