


Anti-parasitic agent 3

Chemical Properties

CAS No. :	2366983-10-8
Formula:	C ₂₁ H ₂₂ FN ₅ O ₃
Molecular Weight:	411.43
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Anti-parasitic agent 3 is effective against drug-resistant parasites.
Targets(IC50)	Others,Parasite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4305 mL	12.1527 mL	24.3055 mL
5 mM	0.4861 mL	2.4305 mL	4.8611 mL
10 mM	0.2431 mL	1.2153 mL	2.4305 mL
50 mM	0.0486 mL	0.2431 mL	0.4861 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481