

## Compound T130440

## Chemical Properties

CAS No. :

Formula: C51H86O28

Molecular Weight: 1147.22

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	2-((6-((2-(((2E,6E,10E)-12-((4-((3,4-dihydroxy-6-methyl-5-((3,4,5-trihydroxy-6-methyloxan-2-yl)oxy)oxan-2-yl)oxy)-3,5-dihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)-3,7,11-trimethyldodeca-2,6,10-trien-1-yl]oxy)-3,5-dihydroxy-6-(hydroxymethyl)oxan-4-yl)oxy]-4,5-dihydroxy-2-methyloxan-3-yl)oxy)-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8717 mL	4.3584 mL	8.7167 mL
5 mM	0.1743 mL	0.8717 mL	1.7433 mL
10 mM	0.0872 mL	0.4358 mL	0.8717 mL
50 mM	0.0174 mL	0.0872 mL	0.1743 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481