

Compound N105-0095

Chemical Properties


CAS No. :

Formula: C₂₅H₂₄O₈

Molecular Weight: 452.459

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N105-0095 is a useful organic compound for research related to life sciences and the catalog number is T131147.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2101 mL	11.0507 mL	22.1014 mL
5 mM	0.442 mL	2.2101 mL	4.4203 mL
10 mM	0.221 mL	1.1051 mL	2.2101 mL
50 mM	0.0442 mL	0.221 mL	0.442 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481