


Pomalidomide-PEG2-OH

Chemical Properties

CAS No. :	2143097-10-1
Formula:	C17H19N3O6
Molecular Weight:	361.35
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Pomalidomide-PEG2-OH enables the synthesis of molecules for targeted protein degradation and PROTAC (proteolysis-targeting chimeras) technology.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7674 mL	13.837 mL	27.674 mL
5 mM	0.5535 mL	2.7674 mL	5.5348 mL
10 mM	0.2767 mL	1.3837 mL	2.7674 mL
50 mM	0.0553 mL	0.2767 mL	0.5535 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481