


Spiroakyraside

Chemical Properties

CAS No. :	139390-87-7
Formula:	C ₅₁ H ₇₈ O ₂₄
Molecular Weight:	1075.161
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Spiroakyraside is a Steroidal from Polygonatum orientale.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9301 mL	4.6505 mL	9.3009 mL
5 mM	0.186 mL	0.9301 mL	1.8602 mL
10 mM	0.093 mL	0.465 mL	0.9301 mL
50 mM	0.0186 mL	0.093 mL	0.186 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Yeşilada E, Houghton PJ. Steroidal saponins from the rhizomes of Polygonatum orientale. *Phytochemistry*. 1991;30(10):3405-9. PubMed PMID: 1367791.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481