


IRE1 $\alpha$  kinase-IN-3

## Chemical Properties

CAS No. :	2416223-41-9
Formula:	C <sub>29</sub> H <sub>32</sub> N <sub>6</sub> O <sub>3</sub> S
Molecular Weight:	544.67
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	IRE1 $\alpha$ kinase-IN-3 is a potent, ATP-competitive inhibitor of IRE1 $\alpha$ (K <sub>i</sub> : 480 nM).
Targets(IC <sub>50</sub> )	Others,IRE1

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.836 mL	9.1799 mL	18.3597 mL
5 mM	0.3672 mL	1.836 mL	3.6719 mL
10 mM	0.1836 mL	0.918 mL	1.836 mL
50 mM	0.0367 mL	0.1836 mL	0.3672 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481