


## O-Desmethyl Mebeverine acid-D5

## Chemical Properties

CAS No. : 1329488-46-1  
 Formula: C<sub>15</sub>H<sub>23</sub>NO<sub>3</sub>  
 Molecular Weight: 270.38  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	O-Desmethyl Mebeverine acid-D5 is the deuterium labeled O-desmethyl Mebeverine acid (T11988).
Targets(IC50)	Drug Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6985 mL	18.4925 mL	36.985 mL
5 mM	0.7397 mL	3.6985 mL	7.397 mL
10 mM	0.3698 mL	1.8492 mL	3.6985 mL
50 mM	0.074 mL	0.3698 mL	0.7397 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481