


Gestodene

Chemical Properties

CAS No. :	60282-87-3
Formula:	C ₂₁ H ₂₆ O ₂
Molecular Weight:	310.43
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Gestodene (WL 70) is a progestogen hormonal contraceptive.
Targets(IC50)	Estrogen/progestogen Receptor, Progesterone Receptor

Solubility Information

Solubility	Ethanol: 8 mg/mL (25.77 mM), Sonication is recommended. DMSO: 58 mg/mL (186.84 mM), Sonication is recommended. H ₂ O: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (6.44 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2213 mL	16.1067 mL	32.2134 mL
5 mM	0.6443 mL	3.2213 mL	6.4427 mL
10 mM	0.3221 mL	1.6107 mL	3.2213 mL
50 mM	0.0644 mL	0.3221 mL	0.6443 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Lidegaard Ø, et al. BMJ. 2011 Oct 25;343:d6423.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481