

## 3-Epi-alpha-amyrin

## Chemical Properties


CAS No. : 5937-48-4

Formula: C<sub>30</sub>H<sub>50</sub>O

Molecular Weight: 426.729

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Epi-alpha-amyrin is a useful organic compound for research related to life sciences. The catalog number is T124591 and the CAS number is 5937-48-4.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3434 mL	11.717 mL	23.434 mL
5 mM	0.4687 mL	2.3434 mL	4.6868 mL
10 mM	0.2343 mL	1.1717 mL	2.3434 mL
50 mM	0.0469 mL	0.2343 mL	0.4687 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481