


2-Mercaptoethyl decanoate

Chemical Properties

CAS No. : 68928-33-6
Formula: C₁₂H₂₄O₂S
Molecular Weight: 232.38
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Mercaptoethyl decanoate is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3033 mL	21.5165 mL	43.033 mL
5 mM	0.8607 mL	4.3033 mL	8.6066 mL
10 mM	0.4303 mL	2.1516 mL	4.3033 mL
50 mM	0.0861 mL	0.4303 mL	0.8607 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481