


S6K-18

Chemical Properties

CAS No. : 1265789-88-5
 Formula: C17H18N4O3S
 Molecular Weight: 358.41
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	S6K-18, a highly selective S6K1 inhibitor, inhibits S6K1 with IC50 of 52nM.
Targets(IC50)	S6 Kinase

Solubility Information

Solubility	DMSO: 3.59 mg/mL (10.02 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7901 mL	13.9505 mL	27.901 mL
5 mM	0.558 mL	2.7901 mL	5.5802 mL
10 mM	0.279 mL	1.3951 mL	2.7901 mL
50 mM	0.0558 mL	0.279 mL	0.558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481