


## Schisphenin E

## Chemical Properties

CAS No. :	1311376-52-9
Formula:	C <sub>24</sub> H <sub>32</sub> O <sub>7</sub>
Molecular Weight:	432.513
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Schisphenin E is a natural product for research related to life sciences. The catalog number is TN6630 and the CAS number is 1311376-52-9.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3121 mL	11.5604 mL	23.1209 mL
5 mM	0.4624 mL	2.3121 mL	4.6242 mL
10 mM	0.2312 mL	1.156 mL	2.3121 mL
50 mM	0.0462 mL	0.2312 mL	0.4624 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481