


## A2A/A1 AR antagonist-1

## Chemical Properties

CAS No. :	2445615-24-5
Formula:	C <sub>15</sub> H <sub>12</sub> N <sub>6</sub> O <sub>2</sub>
Molecular Weight:	308.29
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	A2A/A1 AR antagonist-1 (compound 1a) is a potent dual antagonist with Ki values of 5.58 nM and 24.2 nM for A2A and A1 AR, respectively, demonstrating potential for ischemic stroke research [1].
Targets(IC50)	Others, Adenosine Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2437 mL	16.2185 mL	32.437 mL
5 mM	0.6487 mL	3.2437 mL	6.4874 mL
10 mM	0.3244 mL	1.6218 mL	3.2437 mL
50 mM	0.0649 mL	0.3244 mL	0.6487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481