


SG3-179

Chemical Properties

CAS No. : 1877286-78-6
 Formula: C₂₈H₃₅ClFN₇O₃S
 Molecular Weight: 604.14
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SG3-179 is a BET inhibitor.
Targets(IC50)	Others,FLT,JAK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6552 mL	8.2762 mL	16.5525 mL
5 mM	0.331 mL	1.6552 mL	3.3105 mL
10 mM	0.1655 mL	0.8276 mL	1.6552 mL
50 mM	0.0331 mL	0.1655 mL	0.331 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481