


## Barrinic acid

## Chemical Properties

CAS No. :	38290-00-5
Formula:	C <sub>30</sub> H <sub>46</sub> O <sub>7</sub>
Molecular Weight:	518.691
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Barrinic acid is a useful organic compound for research related to life sciences. The catalog number is T124923 and the CAS number is 38290-00-5.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9279 mL	9.6397 mL	19.2793 mL
5 mM	0.3856 mL	1.9279 mL	3.8559 mL
10 mM	0.1928 mL	0.964 mL	1.9279 mL
50 mM	0.0386 mL	0.1928 mL	0.3856 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481