

GSK-3008348

## Chemical Properties


CAS No. : 1629249-33-7

Formula: C<sub>29</sub>H<sub>37</sub>N<sub>5</sub>O<sub>2</sub>

Molecular Weight: 487.64

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	GSK-3008348 is an antagonist of integrin alpha(v)beta6.
Targets(IC50)	Others,Integrin,TGF-beta/Smad

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0507 mL	10.2535 mL	20.5069 mL
5 mM	0.4101 mL	2.0507 mL	4.1014 mL
10 mM	0.2051 mL	1.0253 mL	2.0507 mL
50 mM	0.041 mL	0.2051 mL	0.4101 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481