

SC 23133

Chemical Properties


CAS No. : 40574-52-5

Formula: C₂₃H₃₀O₃

Molecular Weight: 354.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	SC 23133 is an antimineralocorticoid that is known to induce a marked but reversible inhibition of aldosterone biosynthesis.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.821 mL	14.1052 mL	28.2103 mL
5 mM	0.5642 mL	2.821 mL	5.6421 mL
10 mM	0.2821 mL	1.4105 mL	2.821 mL
50 mM	0.0564 mL	0.2821 mL	0.5642 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481