


(2-Phenylethyl)thio]acetic acid

Chemical Properties

CAS No. :	99186-25-1
Formula:	C10H12O2S
Molecular Weight:	196.266
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2-Phenylethyl)thio]acetic acid ,with CAS No. 99186-25-1, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (2-Phenylethyl)thio]acetic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.095 mL	25.4751 mL	50.9502 mL
5 mM	1.019 mL	5.095 mL	10.190 mL
10 mM	0.5095 mL	2.5475 mL	5.095 mL
50 mM	0.1019 mL	0.5095 mL	1.019 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481