


Pivalanilide

Chemical Properties

CAS No. :	6625-74-7
Formula:	C ₁₁ H ₁₅ NO
Molecular Weight:	177.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pivalanilide is an agent of biochemical.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	-----------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6421 mL	28.2103 mL	56.4207 mL
5 mM	1.1284 mL	5.6421 mL	11.2841 mL
10 mM	0.5642 mL	2.821 mL	5.6421 mL
50 mM	0.1128 mL	0.5642 mL	1.1284 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481