


Azobenzene

Chemical Properties

CAS No. :	103-33-3
Formula:	C ₁₂ H ₁₀ N ₂
Molecular Weight:	182.22
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azobenzene is a phototrigger that can be used to design and synthesize a variety of photoresponsive systems. Azobenzene is genotoxic and induces aggressive sarcomas of the spleen and other abdominal organs.
Targets(IC50)	Others, ROS

Solubility Information

Solubility	DMSO: 150 mg/mL (823.18 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4879 mL	27.4394 mL	54.8787 mL
5 mM	1.0976 mL	5.4879 mL	10.9757 mL
10 mM	0.5488 mL	2.7439 mL	5.4879 mL
50 mM	0.1098 mL	0.5488 mL	1.0976 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wu L, et al. Bioconjug Chem. 2015 Jun 17;26(6):1070-9.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481