


Dodonolide

Chemical Properties

CAS No. :	349534-73-2
Formula:	C ₂₀ H ₂₄ O ₃
Molecular Weight:	312.4
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dodonolide is a natural product of Dodonaea, Sapindaceae. The catalog number is TN3881 and the CAS number is 349534-73-2. Dodonolide can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.201 mL	16.0051 mL	32.0102 mL
5 mM	0.6402 mL	3.201 mL	6.402 mL
10 mM	0.3201 mL	1.6005 mL	3.201 mL
50 mM	0.064 mL	0.3201 mL	0.6402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481