

(2S,3R,4S,5S,6R)-2-[[2-ethoxy-4,5-dihydroxyoxan-3-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol

## Chemical Properties

CAS No. :

Formula: C<sub>13</sub>H<sub>24</sub>O<sub>10</sub>

Molecular Weight: 340.33

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(2S,3R,4S,5S,6R)-2-[[[(2R,3R,4S,5R)-2-ethoxy-4,5-dihydroxyoxan-3-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9383 mL	14.6916 mL	29.3832 mL
5 mM	0.5877 mL	2.9383 mL	5.8766 mL
10 mM	0.2938 mL	1.4692 mL	2.9383 mL
50 mM	0.0588 mL	0.2938 mL	0.5877 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481