


## 5-Acetyltaxachitriene A

## Chemical Properties

CAS No. :	187988-48-3
Formula:	C <sub>34</sub> H <sub>46</sub> O <sub>14</sub>
Molecular Weight:	678.728
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-Acetyltaxachitriene A is a natural product from <i>Taxus mairei</i> .
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4733 mL	7.3667 mL	14.7334 mL
5 mM	0.2947 mL	1.4733 mL	2.9467 mL
10 mM	0.1473 mL	0.7367 mL	1.4733 mL
50 mM	0.0295 mL	0.1473 mL	0.2947 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481