


## Isonemerosin

## Chemical Properties

|                   |   |
|-------------------|---|
| CAS No. :         | 181524-79-8   |
| Formula:          | C <sub>22</sub> H <sub>22</sub> O <sub>7</sub>  |
| Molecular Weight: | 398.41  |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br>Actual storage temperature shall be subject to the COA. |



## Biological Description

|             |  |
|-------------|--|
| Description | Isonemerosin is a natural product for research related to life sciences. The catalog number is TN6592 and the CAS number is 181524-79-8. |
|-------------|--|

## Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 2.510 mL  | 12.5499 mL | 25.0998 mL |
| 5 mM  | 0.502 mL  | 2.510 mL   | 5.020 mL   |
| 10 mM | 0.251 mL  | 1.255 mL   | 2.510 mL   |
| 50 mM | 0.0502 mL | 0.251 mL   | 0.502 mL   |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481