

m-PEG-mal (MW 30000)

**Chemical Properties**

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	m-PEG-mal (MW 30000) is a PEG-based linker for PROTACs, joining two essential ligands crucial for forming PROTAC molecules, and facilitates selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs are composed of two distinct ligands joined by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets the desired protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481