


## Tetrachyrin

## Chemical Properties

CAS No. :	73483-88-2
Formula:	C <sub>20</sub> H <sub>28</sub> O <sub>2</sub>
Molecular Weight:	300.442
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tetrachyrin is a natural product of Wedelia, Asteraceae. The catalog number is TN5128 and the CAS number is 73483-88-2. Tetrachyrin can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3285 mL	16.6423 mL	33.2845 mL
5 mM	0.6657 mL	3.3285 mL	6.6569 mL
10 mM	0.3328 mL	1.6642 mL	3.3285 mL
50 mM	0.0666 mL	0.3328 mL	0.6657 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481