


Y-29794

Chemical Properties

CAS No. : 2554621-98-4
 Formula: C₂₂H₃₂N₂O₅
 Molecular Weight: 404.63
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Y-29794 (Compound 2) is a potent, covalent inhibitor of prolyl oligopeptidase (POP) with a K_i of 0.95 nM, holding research value in neurodegenerative diseases and cancer.
Targets(IC50)	Others,Prolyl Endopeptidase (PREP)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4714 mL	12.357 mL	24.7139 mL
5 mM	0.4943 mL	2.4714 mL	4.9428 mL
10 mM	0.2471 mL	1.2357 mL	2.4714 mL
50 mM	0.0494 mL	0.2471 mL	0.4943 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481