


Izumerogant

Chemical Properties

CAS No. :	2299252-72-3
Formula:	C ₂₂ H ₁₈ ClF ₄ N ₅ O ₂
Molecular Weight:	495.86
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Izumerogant is an inverse agonist of the retinoid-related orphan receptor-gamma (ROR γ).
Targets(IC50)	Others,ROR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0167 mL	10.0835 mL	20.167 mL
5 mM	0.4033 mL	2.0167 mL	4.0334 mL
10 mM	0.2017 mL	1.0083 mL	2.0167 mL
50 mM	0.0403 mL	0.2017 mL	0.4033 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481