

CytoCalcein Violet 500

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CytoCalcein Violet 500 is a fluorogenic dye used to assess cell viability. It is a cell-permeable dye that, upon entering cells, is cleaved by intracellular esterases. CytoCalcein Violet 500 is retained by living cells and displays excitation/emission maxima of 405/500 nm, respectively.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481