

## Senampeline B

## Chemical Properties


CAS No. :

Formula: C<sub>25</sub>H<sub>31</sub>N<sub>0</sub>S<sub>0</sub>

Molecular Weight: 473.522

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Senampeline B is a useful organic compound for research related to life sciences and the catalog number is T124518.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1118 mL	10.5592 mL	21.1184 mL
5 mM	0.4224 mL	2.1118 mL	4.2237 mL
10 mM	0.2112 mL	1.0559 mL	2.1118 mL
50 mM	0.0422 mL	0.2112 mL	0.4224 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481