


Apazone dihydrate

Chemical Properties

CAS No. :	22304-30-9
Formula:	C ₁₆ H ₂₄ N ₄ O ₄
Molecular Weight:	336.39
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Apazone dihydrate is a nonsteroidal anti-inflammatory agent.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9727 mL	14.8637 mL	29.7274 mL
5 mM	0.5945 mL	2.9727 mL	5.9455 mL
10 mM	0.2973 mL	1.4864 mL	2.9727 mL
50 mM	0.0595 mL	0.2973 mL	0.5945 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Kim HS, Austin J, Hage DS. Screening major binding sites on human serum albumin by affinity capillary electrophoresis. *Methods Mol Biol.* 2004;276:169-87. PubMed PMID: 15163858.
- García Rodríguez LA, Hernández-Díaz S. Relative risk of upper gastrointestinal complications among users of acetaminophen and nonsteroidal anti-inflammatory drugs. *Epidemiology.* 2001 Sep;12(5):570-6. PubMed PMID: 11505178.
- Géher P, Gömör B. Repeated cyclosporine therapy of peripheral arthritis associated with ankylosing spondylitis. *Med Sci Monit.* 2001 Jan-Feb;7(1):105-7. PubMed PMID: 11208503.
- Saenz A, Ausejo M, Shea B, Wells G, Welch V, Tugwell P. Pharmacotherapy for Behcet's syndrome. *Cochrane Database Syst Rev.* 2000;(2):CD001084. Review. PubMed PMID: 10796413.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481