

Propargyl-O-C1-amido-PEG4-C2-NHS ester

Chemical Properties

CAS No. : 2101206-92-0


Formula: C₂₀H₃₀N₂O₁₀

Molecular Weight: 458.46

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|--|
| Description | Propargyl-O-C1-amido-PEG4-C2-NHS ester is a non-cleavable 4-unit PEG linker used in antibody-drug conjugation (ADC) to link antibodies with drugs. |
| Targets(IC50) | ADC Linker,PROTAC Linker |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.1812 mL | 10.9061 mL | 21.8122 mL |
| 5 mM | 0.4362 mL | 2.1812 mL | 4.3624 mL |
| 10 mM | 0.2181 mL | 1.0906 mL | 2.1812 mL |
| 50 mM | 0.0436 mL | 0.2181 mL | 0.4362 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481