


Ginsenoside Re

Chemical Properties

CAS No. :	52286-59-6
Formula:	C ₄₈ H ₈₂ O ₁₈
Molecular Weight:	947.15
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ginsenoside Re (Ginsenoside B2) may have properties that inhibit or prevent the growth of tumors.
Targets(IC50)	Beta Amyloid,NF-κB,Endogenous Metabolite,JNK

Solubility Information

Solubility	DMSO: 250 mg/mL (263.95 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (2.11 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0558 mL	5.279 mL	10.558 mL
5 mM	0.2112 mL	1.0558 mL	2.1116 mL
10 mM	0.1056 mL	0.5279 mL	1.0558 mL
50 mM	0.0211 mL	0.1056 mL	0.2112 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Nam Y, et al. J Appl Toxicol. 2015 Aug;35(8):927-44.

Chen Z, Ni R, Hu Y, et al. A natural protopanaxatriol from Panax notoginseng enhances osteosarcoma sensitivity to ferroptosis via ASCL4 upregulation. Journal of Functional Foods. 2024, 122: 106488.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481