


## Benzyl-PEG7-alcohol

## Chemical Properties

CAS No. : 423763-19-3

Formula: C<sub>21</sub>H<sub>36</sub>O<sub>8</sub>

Molecular Weight: 416.51



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Benzyl-PEG7-alcohol is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binds to an E3 ubiquitin ligase and the other to the target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4009 mL	12.0045 mL	24.009 mL
5 mM	0.4802 mL	2.4009 mL	4.8018 mL
10 mM	0.2401 mL	1.2005 mL	2.4009 mL
50 mM	0.048 mL	0.2401 mL	0.4802 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481