


## Ehp-inhibitor-2

## Chemical Properties

CAS No. :	861249-77-6
Formula:	C17H13N5O
Molecular Weight:	303.32
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ehp-inhibitor-2 is an Eph family tyrosine kinase inhibitor targeting Eph receptors.
Targets(IC50)	Ephrin Receptor

## Solubility Information

Solubility	DMSO: 27.5 mg/mL (90.66 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	----------------------------------------------------------------------------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2968 mL	16.4842 mL	32.9685 mL
5 mM	0.6594 mL	3.2968 mL	6.5937 mL
10 mM	0.3297 mL	1.6484 mL	3.2968 mL
50 mM	0.0659 mL	0.3297 mL	0.6594 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481