


## 11(R)-HEPE

## Chemical Properties

CAS No. :	109430-11-7
Formula:	C <sub>20</sub> H <sub>30</sub> O <sub>3</sub>
Molecular Weight:	318.457
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	11(R)-HEPE is produced by the oxidation of EPA by 11(R)-LO. This enzymatic activity and the resulting 11(R)-hydroxy acid have been isolated from the sea urchin <i>S. purpuratus</i> .
Targets(IC50)	Others

## Solubility Information

Solubility	0.1 M Na <sub>2</sub> CO <sub>3</sub> : 2 mg/mL (6.28 mM), Sonication is recommended. DMF: Miscible Ethanol: Miscible DMSO: Miscible PBS (pH 7.2): 0.8 mg/mL (2.51 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1401 mL	15.7006 mL	31.4011 mL
5 mM	0.628 mL	3.1401 mL	6.2802 mL
10 mM	0.314 mL	1.5701 mL	3.1401 mL
50 mM	0.0628 mL	0.314 mL	0.628 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481