

Ethyl 2-(pyridin-3-yl)acetate

Chemical Properties


CAS No. : 39931-77-6

Formula: C₉H₁₁NO₂

Molecular Weight: 165.19

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethyl 2-(pyridin-3-yl)acetate, with CAS No. 39931-77-6, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Ethyl 2-(pyridin-3-yl)acetate provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0536 mL	30.2682 mL	60.5364 mL
5 mM	1.2107 mL	6.0536 mL	12.1073 mL
10 mM	0.6054 mL	3.0268 mL	6.0536 mL
50 mM	0.1211 mL	0.6054 mL	1.2107 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481