


## 6-Methoxynaringenin

## Chemical Properties

CAS No. :	94942-49-1
Formula:	C <sub>16</sub> H <sub>14</sub> O <sub>6</sub>
Molecular Weight:	302.28
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	6-Methoxynaringenin is an aerial part of <i>Scutellaria barbata</i> D. Don (Lamiaceae) and inhibits NO production efficiently with an IC <sub>50</sub> value of 25.8 μM.
Targets(IC <sub>50</sub> )	NO Synthase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3082 mL	16.541 mL	33.0819 mL
5 mM	0.6616 mL	3.3082 mL	6.6164 mL
10 mM	0.3308 mL	1.6541 mL	3.3082 mL
50 mM	0.0662 mL	0.3308 mL	0.6616 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Lee SR, et al. Constituents from *Scutellaria barbata* Inhibiting Nitric Oxide Production in LPS-Stimulated Microglial Cells. *Chem Biodivers.* 2017 Nov;14(11).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481