


## 2'-C-Methyl-6-O-methylinosine

## Chemical Properties

CAS No. :	565450-78-4
Formula:	C <sub>12</sub> H <sub>16</sub> N <sub>4</sub> O <sub>5</sub>
Molecular Weight:	296.28
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2'-C-Methyl-6-O-methylinosine (2'-C-Methyl-6-O-methylinosine) is a hypoxanthine analog exhibiting diverse physical activities for researching various tumor diseases.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Xanthine Oxidase

## Solubility Information

Solubility	DMSO: 60 mg/mL (202.51 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3752 mL	16.8759 mL	33.7519 mL
5 mM	0.675 mL	3.3752 mL	6.7504 mL
10 mM	0.3375 mL	1.6876 mL	3.3752 mL
50 mM	0.0675 mL	0.3375 mL	0.675 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481