


Lu AA 47070

Chemical Properties

CAS No. : 913842-25-8
 Formula: C₁₇H₂₀F₂N₃O₆P₅
 Molecular Weight: 463.39
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	adenosine A2A receptor antagonist
Targets(IC50)	Others,Adenosine Receptor

Solubility Information

Solubility	DMSO: <46.34 mg/mL,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.158 mL	10.790 mL	21.5801 mL
5 mM	0.4316 mL	2.158 mL	4.316 mL
10 mM	0.2158 mL	1.079 mL	2.158 mL
50 mM	0.0432 mL	0.2158 mL	0.4316 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481