

## Acid-PEG1-C2-Boc

## Chemical Properties

CAS No. : 2086689-08-7


Formula: C<sub>10</sub>H<sub>18</sub>O<sub>5</sub>

Molecular Weight: 218.25

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Acid-PEG1-C2-Boc is a PEG-based linker for PROTACs that joins two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands joined by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets the protein of interest. They leverage the intracellular ubiquitin-proteasome system for the selective degradation of target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5819 mL	22.9095 mL	45.819 mL
5 mM	0.9164 mL	4.5819 mL	9.1638 mL
10 mM	0.4582 mL	2.291 mL	4.5819 mL
50 mM	0.0916 mL	0.4582 mL	0.9164 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481