

## 1-(2,3-dihydroxypropoxy)-1-oxononan-3-yl 3-hydroxynonanoate

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>40</sub>O<sub>7</sub>

Molecular Weight: 404.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	1-(2,3-dihydroxypropoxy)-1-oxononan-3-yl 3-hydroxynonanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4719 mL	12.3597 mL	24.7194 mL
5 mM	0.4944 mL	2.4719 mL	4.9439 mL
10 mM	0.2472 mL	1.236 mL	2.4719 mL
50 mM	0.0494 mL	0.2472 mL	0.4944 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481