


16-Dehydroprogesterone

Chemical Properties

CAS No. :	1096-38-4
Formula:	C ₂₁ H ₂₈ O ₂
Molecular Weight:	312.45
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	16-Dehydroprogesterone is an endogenous Metabolite. 16-Dehydroprogesterone is a steroidal progestin.
Targets(IC50)	Endogenous Metabolite

Solubility Information

Solubility	DMSO: 14.71 mg/mL (47.08 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (6.4 mM), Sonication is recommended. 10% DMSO+90% Saline: 1.47 mg/mL (4.7 mM), Solution. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2005 mL	16.0026 mL	32.0051 mL
5 mM	0.6401 mL	3.2005 mL	6.401 mL
10 mM	0.3201 mL	1.6003 mL	3.2005 mL
50 mM	0.064 mL	0.3201 mL	0.6401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Glass TL, et al. Stimulation of 16-dehydroprogesterone and progesterone reductases of

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481