


EP 157

Chemical Properties

CAS No. : 101910-65-0
 Formula: C₂₉H₃₃NO₃
 Molecular Weight: 443.587
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|---------------------------------|
| Description | EP 157 is a bioactive chemical. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.2543 mL | 11.2717 mL | 22.5433 mL |
| 5 mM | 0.4509 mL | 2.2543 mL | 4.5087 mL |
| 10 mM | 0.2254 mL | 1.1272 mL | 2.2543 mL |
| 50 mM | 0.0451 mL | 0.2254 mL | 0.4509 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481